**EN SARS-CoV-2 & Influenza A/B & RSV Antigen Kit  
(Colloidal Gold)**

**Catalog No.** CG127001, CG127005, CG127025

**1. Intended Use**

This SARS-CoV-2 & Influenza A/B & RSV Antigen Kit (Colloidal Gold) is a lateral flow rapid chromatographic immunoassay intended for the in vitro qualitative detection and differentiation of antigens from SARS-CoV-2, influenza A, influenza B and Respiratory Syncytial Virus (RSV) in nasal, nasopharyngeal, and oropharyngeal (throat) swab specimens from individuals suspected of respiratory viral infections.

Results are for the identification of antigens of SARS-CoV-2, influenza A, influenza B and RSV, but does not differentiate between SARS-CoV and SARS-CoV-2 viruses and is not intended to detect influenza C antigens. These viral antigens are generally detectable in upper respiratory specimens during the acute phase of infection. Positive results indicate the presence of viral antigens, but correlation with patient history and other diagnostic information is necessary to determine the infection status. Positive results do not rule out bacterial infection or co-infection with other viruses. The agent detected may not be the definite cause of disease.

Negative results should be treated as presumptive, do not rule out either influenza, SARS-CoV-2 or RSV infection, and should not be used as the sole basis for treatment or patient management decisions, including infection control decisions. It's recommended to confirm negative SARS-CoV-2 results with a molecular assay, and negative influenza A/B and RSV results with viral culture or a molecular assay, if necessary, for patient management.

**2. Test Principle**

Colloidal gold based immunochromatography is applied. The test is a combination of SARS-CoV-2 antigen test, influenza A/B antigen test and RSV antigen test. Briefly, the SARS-CoV-2 or influenza A/B or RSV antigens in the sample bind with the colloidal gold (detector) labeled anti-SARS-CoV-2 / anti-influenza A/B / anti-RSV antibodies on the conjugate pad to form the immunocomplexes. When the complexes migrate to the test line (coated with anti-SARS-CoV-2 or influenza A/B or RSV antibodies), the complexes that contain the SARS-CoV-2 or influenza A/B or RSV antigen will be captured. Similarly, the colloidal gold labeled chicken IgY antibody will be captured in the control lines (coated with goat anti-chicken IgY). The complexes containing SARS-CoV-2 or influenza A/B or RSV antigen will produce a colored line in the specific test line region, indicating the sample is positive for SARS-CoV-2 or influenza A/B or RSV antigens. Similarly, a colored line will also appear in each of the control line region indicating the adequate volume of specimen has been added and the membrane wicking has occurred.

**3. Materials Provided**

Component	1T CG127001	5T CG127005	25T CG127025
Test Cartridge	1	5	25
Sampling Swab	1	5	25
Extraction Tube	1	5	25
Instructions For Use	1	1	1

**4. Material Required but Not Supplied**

4.1 Timer

4.2 Any necessary personal protective equipment

**5. Storage and Stability**

Store the kit at 2 – 30°C in a dry place and avoid direct sunlight. Do not freeze. The unopened cartridges are stable until the expiry date printed on the labels. It's recommended to open the test cartridge pouch just before use. Once the pouch is opened, the test should be started within 60 minutes.

**6. Sample Collection**

6.1 Nasal, nasopharyngeal, and oropharyngeal (throat) swab samples are acceptable for testing with this kit.

6.2 **Proper sample collection and handling is critical to the performance of the kit.** When collecting the swab sample, use the swab supplied in the kit. Do NOT touch the swab tip when handling the swab sample.

(a). Nasal swab collection: while gently rotating, insert the swab into one nostril. The swab tip should be inserted up to 2.5 cm (1 inch) from the edge of the nostril. Rotate the swab 5 times against the nasal wall. Remove and repeat the sampling process using the same swab for the other nostril.

(b). Nasopharyngeal swab collection: insert swab through the nares perpendicular to the nose (face) until resistance is encountered and the

fingers touch the nose. Leave the swab in place for 15 – 30 seconds. Rotate the swab 3 times and remove it from the nasopharynx.

(c). Oropharyngeal (throat) swab collection: insert swab into the posterior pharynx and tonsillar areas. Rub swab over both tonsillar pillars and posterior oropharynx and avoid touching the tongue, teeth, and gums.

6.3 Samples should be tested immediately after collection.

**7. Test Procedure**

7.1 Allow the kit to come to room temperature (15 – 30°C) before testing.

7.2 Remove the Test Cartridge from its package. Place the cartridge on a clean, flat and dry surface. Label the cartridge with patient ID.

7.3 Remove the Dropper Lid from the Extraction Tube. Remove the cap and the inner plug (stopper) from the Dropper Lid.

7.4 Place and soak the patient swab into the Extraction Tube. Break the swab stick at the breaking point to leave the lower half of the swab in the Extraction Tube.

7.5 Attach the Dropper Lid firmly onto the top of the Extraction Tube. Squeeze the swab from outside the tube 10 times. Leave the swab in the tube for 1 minute.

7.6 Invert the Extraction Tube, and then vertically add **2 drops** into each sample well of the cartridge.

7.7 Start the timer: leave the cartridge at room temperature for **15 minutes**.

7.8 After the 15-minute incubation, read the results according to section 8. **Interpretation of Results.** Do not interpret the results after 20 minutes (from addition of the sample).

Note: Make sure there is sufficient light when reading and interpreting the results. Results read before 15 minutes or after 20 minutes may be false negative, false positive, or invalid.

**8. Interpretation of Results**

For each of the test lane, i.e., SARS-CoV-2, Influenza A/B and RSV tests, the Control (C) line must appear for the test to be valid regardless of the appearance of the test line(s). Failure to observe the control line indicates the results are not reliable. When this occurs, check the operation procedure carefully, and test again with a new sample. If the problem recurs, contact your local distributor.

**For SARS-CoV-2, read the results on the first window.**

**For Influenza A/B, read the results on the second window.**

**For RSV, read the results on the third window.**

Results	SARS-CoV-2 Result Window Lane 1	Influenza A/B Result Window Lane 2	RSV Result Window Lane 3
Positive <sup>1)2)</sup>	 SARS-CoV-2 Positive	 Influenza A Positive   Influenza B Positive   Influenza A & B Positive	 RSV Positive
Negative <sup>3)</sup>	 SARS-CoV-2 Negative	 Influenza Negative	 RSV Negative
Invalid <sup>4)</sup>	 Invalid SARS-CoV-2 Test	 Invalid Influenza Test	 Invalid RSV Test

1) The test lines (A, B or T line region) may vary in intensity depending on the amount of antigen detected. **Any visible (even light or faint) T line(s) should be interpreted as positive.**

2) Co-infection with influenza A, B and/or SARS-CoV-2 and/or RSV is rare. If the results are positive for more than one antigen, e.g., influenza A, B and/or SARS-CoV-2, it's recommended to retest with a new sample.

3) Negative results should be treated as presumptive. Negative results do not rule out influenza or SARS-CoV-2 or RSV infection, and should not be used as the sole basis for treatment or patient management decisions.

4) Insufficient sample volume or incorrect operation are the most likely reasons for invalid results. Review the procedure and repeat the test with a new sample.

**9. Limitations**

- 9.1 The test is only for qualitative detection of influenza A, influenza B, SARS-CoV-2 and RSV antigens.
- 9.2 This test detects both viable and non-viable influenza A, influenza B, SARS-CoV-2 and RSV. Test performance depends on the amount of virus (antigen) in the sample and may or may not correlate with viral culture or molecular assay performed on the same sample.
- 9.3 A negative test result may occur if the level of antigen in a sample is below the detection limit of the test, if the sample was collected or handled improperly or if the virus has amino acid changes in the target epitope region recognized by the antibodies utilized in the test.
- 9.4 Failure to follow the Test Procedure may adversely affect test performance and/or invalidate the test result.
- 9.5 Test results obtained with the test must be used in conjunction with other clinical history, findings, and epidemiological data, etc.
- 9.6 Positive test results do not rule out co-infections with other pathogens, differentiate between SARS-CoV and SARS-CoV-2 or identify specific influenza subtypes.
- 9.7 Negative results should be treated as presumptive, do not rule out either influenza or SARS-CoV-2 or RSV infection, and should not be used as the sole basis for treatment or patient management decisions, including infection control decisions. Negative results must be combined with clinical observations, patient history and epidemiological information, and confirmed with a molecular assay, if necessary, for patient management.

**10. Performance**

**10.1 Clinical Performance**

The clinical performance of SARS-CoV-2 & Influenza A/B & RSV Antigen Kit (Colloidal Gold) was determined by comparing to RT-PCR tests for SARS-CoV-2, influenza A and B and RSV using clinical samples.

Test	Sample Type	N	Sensitivity	Specificity
SARS-CoV-2	Nasal swab	200	93.00%	100.00%
	Nasopharyngeal swab	268	94.05%	100.00%
	Oropharyngeal	200	91.00%	100.00%
	Combined	668	92.93%	100.00%
Influenza A	Nasal swab	159	85.00%	99.28%
	Nasopharyngeal swab	86	95.24%	100.00%
	Oropharyngeal	70	95.00%	100.00%
	Combined	315	91.80%	99.61%
Influenza B	Nasal swab	159	95.00%	100.00%
	Nasopharyngeal swab	86	100.00%	100.00%
	Oropharyngeal	71	85.00%	100.00%
	Combined	316	93.44%	100.00%
RSV	Nasal swab	128	92.31%	100.00%
	Nasopharyngeal swab	120	100.00%	100.00%
	Oropharyngeal	75	100.00%	100.00%
	Combined	323	96.97%	100.00%

**10.2 Limit of Detection (LoD)**

The kit was confirmed to detect 49 TCID<sub>50</sub>/mL of SARS-CoV-2, 4 x 10<sup>4</sup> TCID<sub>50</sub>/mL of seasonal influenza A(H3N2), 2 x 10<sup>4</sup> TCID<sub>50</sub>/mL of seasonal influenza A(H1N1), 1.96 x 10<sup>4</sup> TCID<sub>50</sub>/mL of influenza A(H1N1) pdm09, 2.63 x 10<sup>3</sup> TCID<sub>50</sub>/mL of influenza B (Yamagata), 5 x 10<sup>3</sup> TCID<sub>50</sub>/mL of influenza B (Victoria), and 1.6 x 10<sup>4</sup> TCID<sub>50</sub>/mL of RSV.

**10.3 Cross Reactivity and microbial interference**

The cross-reactivity and microbial interference study evaluated whether microorganisms possibly contained in clinical samples interfere with the kit which may lead to false positive or false negative results. Each microorganism was tested in triplicate in the presence (3 x LoD) or absence of the analytes. No cross-reactivity or interference with the microorganisms listed in the table below was found.

No.	Microorganism	Final Test Concentration
1	HCoV-OC43	2 x 10 <sup>6</sup> TCID <sub>50</sub> /mL
2	HCoV-229E	2 x 10 <sup>6</sup> TCID <sub>50</sub> /mL

3	HCoV-NL63	2 x 10 <sup>6</sup> TCID <sub>50</sub> /mL
4	Rotavirus	2 x 10 <sup>6</sup> TCID <sub>50</sub> /mL
5	MERS	1 x 10 <sup>6</sup> TCID <sub>50</sub> /mL
6	Adenovirus	2 x 10 <sup>6</sup> TCID <sub>50</sub> /mL
7	Norovirus	2 x 10 <sup>6</sup> TCID <sub>50</sub> /mL
8	Mycoplasma pneumonia	1.5 x 10 <sup>6</sup> cfu/mL

**10.4 Interference**

The following interfering substances have no impact on SARS-CoV-2 & Influenza A/B & RSV Antigen Kit (Colloidal Gold).

No.	Interfering substance	Final Test Concentration
1	Phenylephrine	15% v/v
2	Oxymetazoline	15% v/v
3	Sodium chloride	5 mg/mL
4	Beclomethasone	5 ng/mL
5	Dexamethasone	0.5 µg/mL
6	Flunisolide	0.5 µg/mL
7	Triamcinolone acetonide	1 ng/mL
8	Budesonide	2.5 ng/mL
9	Mometasone	1 ng/mL
10	Fluticasone	2 ng/mL

**11. Warnings and Precautions**

- 11.1 For in vitro diagnostic use only.
- 11.2 If the package has been damaged, the label cannot be seen clearly or if the kit has expired, do not use the kit.
- 11.3 Read the instructions carefully before starting the test. To obtain accurate results, the instructions must be followed.
- 11.4 Extraction Buffer of different lots are not interchangeable. The results may not be reliable if reagents from different lots are mixed or used together.
- 11.5 The test cartridge is for single test and cannot be reused. Do not use expired cartridges.
- 11.6 Do not eat the desiccant.
- 11.7 Avoid exposure of the skin, eyes, mouth, mucous membranes to the Extraction Buffer. If the contact occurs, flush with large amounts of water.
- 11.8 Use of gloves during sampling and testing is recommended.
- 11.9 The samples, used reagents and consumables are potentially infectious waste and should be disposed of in accordance with national and local regulations.

**12. Symbols**

Consult instructions for use	LOT Lot number	Do not reuse
Use-by date	In vitro diagnostic medical device	Temperature limit 2 - 30°C
Manufacturer	Tests per kit	Avoid sunshine
Date of manufacture	REF Catalog number	Authorized representative in the European Community
This product complies with the requirements of the European Directive 98/79/EC		

**Goldsite Diagnostics Inc.**

No. 103C, 503C & 504D, Technology Building & No. 3A & 4A, Technology Building Annex, Zhaoshang Sub-District, Nanshan District, Shenzhen, China, 518067

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